

Unleash the power of conversation

Signia Integrated Xperience

with RealTime Conversation Enhancement



signia

Insio IX
Silk Charge&Go IX
Pure Charge&Go IX
Styletto IX





95%

of wearers have better speech reception with RealTime Conversation Enhancement than without.¹

The science of being part of the conversation

With RealTime Conversation Enhancement,² your hearing aids can now enhance and follow multiple conversation partners simultaneously, even as they move or you turn your head. Available in select models, this technology works in three stages:



Analyze: Evaluates 192,000 data points every second to build a picture of your listening situation and to understand your environment and conversation needs.



Augment: Boosts the voices of all the people you are talking to, helping you hear them clearly while keeping you immersed in your surroundings.



Adapt: Updates 1,000 times a second to your conversation environment allowing it to seamlessly adjust in real-time as the conversation moves and shifts, ensuring you are always part of the action.

¹ Visit signia-pro.com/en-us/claims-data for more details.

² RealTime Conversation Enhancement not available in all styles.




Finally, a hearing aid designed for real conversations

Every conversation is a powerful tool to connect, engage, and bring us closer to realizing our potential. Real conversations do not stand still – they are vibrant, unpredictable, and full of surprises.

Signia Integrated Xperience is built from the ground up with conversation in mind. It is powered by revolutionary technology that helps you get involved and keep up with the conversation, no matter how busy it gets.

Contribute to the conversation with ease, and never miss a beat. With Signia Integrated Xperience, you are not just part of the conversation. You are at its heart.

NEW



Insio IX

Our most discreet custom hearing aid is now reengineered to masterfully reduce background noise and keep you in the conversation. Embrace an almost invisible hearing aid, tailored just for you.

- Binaural OneMic Directionality provides focused conversation, even in noisy environments
- Speech Stabilizer for clearer conversations in dynamic noisy environments
- Reengineered for comfort and our smallest custom device ever
- Custom-molded for better retention, increased comfort and reduced sound leakage reducing feedback compared to an instant-fit device
- Floating mic for additional flexibility in design, allowing a deeper, more secure fitting
- IIC and CIC styles
- 3 color options



Available Faceplate Colors:



Black



Dark
Brown



Mocha



Silk Charge&Go IX

The confidence you get from being able to hear conversations clearly makes it easy to feel more involved. And that feeling grows when your hearing aid is practically invisible and can be fitted in an instant.

Comfortable and discreet, this tiny hearing aid provides the convenience of on-the-go rechargeability. It is packed with advanced sound features to support you in any environment including its unique Binaural OneMic Directionality, which offers unmatched conversation support in noisy environments for a hearing aid of its size.

- Rechargeable with up to 28 hours of battery life
- Ultra-discreet, instant fit
- Speech Stabilizer for clearer conversations in dynamic noisy environments
- Portable charger with 4 days of additional charges



**Available
Faceplate
Colors:**



Black



Mocha



Pure Charge&Go IX

Pure Charge&Go IX hearing aids are the easiest way to enjoy group conversations again. Small, discreet, and conveniently rechargeable for all-day listening, with Bluetooth connectivity to smartphones. Optional Pure Charge&Go T IX offers t-coil and extended battery life.

- Processes speech separately from the surrounding sounds for comfort and clarity
- Up to 39 hours of battery life
- Small and sleek design
- Bluetooth streaming to iPhone® and Android™
- 10 color options
- Own Voice Processing 2.0 to make your voice sound more natural
- Get help when you need it, where you need it with Signia Assistant real-time adjustments directly in the Signia app



Available Colors:



Black



Graphite



Dark
Champagne



Silver



Pearl
White



Fine
Gold



Deep
Brown




Sandy
Brown



Rose
Gold



Beige



Styletto IX

Sleek, sophisticated, and available in various distinctive colors, Styletto hearing aids blend superior IX technology with style. Enjoy great hearing and a great look — and the added convenience of recharging on the go.

- Processes speech separately from the surrounding sounds for comfort and clarity
- Refined ultra-slim design
- Rechargeable for a full day: Styletto IX charger provides up to five full charges before needing to be plugged in
- CROS compatibility
- 7 color options
- Own Voice Processing 2.0 to make your voice sound more natural
- Get help when you need it, where you need it with Signia Assistant real-time adjustments directly in the Signia app



Available Colors:



Black/
Black Gloss



Black/
Graphite



Black/
Silver



Cosmic Blue/
Rose Gold



Snow White/
Rose Gold



Snow White/
Silver



White/
White Gloss

About Signia

Signia is one of the leading hearing aid brands in the world. Our goal is to enhance human performance through iconic innovations that shape the prescription audio category. We refuse to accept hearing loss as a limitation. It's not about correcting a loss, it's about gaining an edge. Not just performing, but performing brilliantly. Not just hearing again, but hearing better than ever. Hearing every word, beat, and breath. With Signia, you can Be Brilliant.

The information in this document contains general descriptions of the technical options available, which do not always have to be present in individual cases and are subject to change without prior notice.

The Bluetooth® word mark and logos are owned by the Bluetooth SIG, Inc. and any use of such marks by WSAUD A/S is under license. Other trademarks and trade names are those of their respective owners. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. Android, Google Play and the Google Play logo are trademarks of Google Inc.



© 02.2024 Signia is a registered trademark of WSAUD A/S SIV-18461

Find us on:



signia.net



facebook.com/signiahearing



youtube.com/signiahearing



instagram.com/signiahearing

**Be
Brilliant™**